Ultra-Ma











storage case				S
Product name	Ultrasonic device for the head	Ultrasonic device for the body	Ultrasonic device for the head and body	Ultrasonic device for the jaw joint
Model	UM-123 / UM-123X	UM-124 ① / ② / ③	Bp-21 ① / ②	Hp-20
Power supply : AC 100-240 V	0.2 A or less	0.5 A or less	0.5 A or less	0.5 A or less
Output frequency	30 kHz ±10%	30 kHz ±20%	30 kHz ±20%	30 kHz ±20%
Maximum ultrasonic output intensity	0-10 mW or less	0-100 mW or less	0-100 mW or less	0-100 mW or less
Music intensity	0~2 mW	0~2 mW	0~2 mW	0~2 mW
Vibration	-	0.5 W (special version)	0.5 W (special version)	-
Timer	Initial setting: 20 min (standard use: 20 min)	Initial setting: 20 min (standard use: 20 min)	Initial setting: 20 min (standard use: 20 min)	Initial setting: 20 min (standard use: 20 min)
Mozart's music in storage	Serenade / Divertimento / German Dance	Serenade / Divertimento / German Dance	Serenade / Divertimento / German Dance	Serenade / Divertimento / German Dance
Port for external audio source	3.5 mm jack	-	3.5 mm jack	3.5 mm jack
Control unit body dimensions	11×6×19 cm (WHD)	32×7×23 cm (WHD)	11×6×19 cm (WHD)	11×6×19 cm (WHD)
Storage case dimensions/ total weight	43×13×28 cm (WHD) 2.5 kg	46×15×33 cm (WHD) 5.4 kg / 5.8 kg / 5.8 kg	43×13×28 cm (WHD) 3.1 kg	43×13×28 cm (WHD) 3.1 kg
Standard price (tax included)	UM-123 ··· ¥682,000 UM-123X ··· ¥880,000	① ··· ¥1,980,000 ② ··· ¥3,080,000 ③ ··· ¥3,190,000	① ··· ¥2,200,000 ② ··· ¥2,300,000	¥1,980,000

- ③ ··· ¥3,190,000 • As of January 2023, prices currently reflect the 10% consumption tax, but in the event of a change in the tax, prices will also change.
- The specifications and design of the products listed may change without prior notice as improvements are made.
- The actual products may differ slightly from the pictures included due to limitation of the printing process.
- The rights to all information concerning these products belong to Kamiyama, Ltd.

Ultra-Ma and Mu-Ma are registered trademarks of Kamiyama, Ltd.

Kamiyama

⟨ Developer/manufacturer ⟩

Kamiyama, Ltd.

2-264 Kamiyamacho, Funabashi, Chiba 273-0046, Japan

TEL·FAX: +81-47-438-0242

- •Please see our website for more information about our products and the posted paper.
- ·In some cases, a password may be necessary. Please contact us.



https://worldbrain.jp

⟨ General distributor ⟩

BP-TECH Inc.

Japan Selfcare Institute®

Head office: 2FL, 2-10-29 Kanda Jimbocho, Chiyoda-ku, Tokyo 101-0051, Japan

You can see and try all of our products at the showroom below. You can also discuss purchases or rentals there.

Brain & Body Store® (weekdays 10:00-17:00)

2-9-13 Azamino, Aoba-ku, Yokohama, Kanagawa 225-0011, Japan

TEL:+81-45-905-3365 FAX:+81-45-905-3369





https://www.bp-tech.jp/ultra-ma (online shop)



https://www.japan-selfcare-i.jp/ultra (online shop and rental))

* Some products are also available on Amazon.



Health maintenance through ultrasound.



An ultrasound approach for your head and body (Ultra-Ma)

Dolphins' ultrasonic waves and Mozart

The same 30 kilohertz ultrasonic waves that dolphins give off grow to a 5 cm wavelength inside the human body, making for a soft approach. What's more, you can listen to Mozart's music through a bone conduction speaker as you relax at your leisure.

Ultra-Ma Ultrasonic device for the head

UM-123 / UM-123X



This is headband-type ultrasonic device.

- You can adjust the strength of the ultrasonic waves and the volume of the music.
- It contains 3 tunes by Mozart.
- You can connect to the external port to listen to any music you like.

Set contents

Control unit

Headband-type Ultrasonic device





Navy blue

UM-123 ··· ¥682,000 (tax included)

〈 Specifications 〉 Ultrasonic transducers are attached in 2 locations on the front portion of the headband.

UM-123X ··· ¥880,000 (tax included)

⟨ Specifications ⟩ Ultrasonic transducers are attached in 2 locations on the front portion of the headband as well as 1 each on the left and right sides.

Ultra-Ma Ultrasonic device for the body UM-124



These are 3 types of bodyband-type ultrasonic device. Use them on areas other than the head.

- The ultrasonic device for the body stimulates with ultrasonic transducers with about 10 times the intensity of the ultrasonic device for the head.
- The ultrasonic device for the body sends strong and weak waves through the body about 6,000 times/second, transmitting micro movements to your anatomy.
- Use the body band by wrapping it directly against the skin, or placing it against the area that concerns you.
- With the ultrasonic device for the body, you can choose from music oscillation (A model) or vibration (B model).
- You can adjust the strength of the ultrasonic waves and the volume of the music.
- It contains 3 tunes by Mozart.



* The bodyband is not worn over clothing. but rather directly against the skin.

Set contents

Control unit



Ultrasonic device Ultrasonic device Music oscillator Vibration





〈 A model 〉 〈 B model 〉

Bodyband



The double band sets include 2 bodybands.

UM-124 ① ··· ¥1,980,000 (tax included) 〈 Contents 〉 1 bodyband-type A model

UM-124 ② ··· ¥3,080,000 (tax included) 〈 Contents 〉 [double bands] 2 bodyband-type A models

UM-124 ③ ··· ¥3,190,000 (tax included) 〈 Contents 〉 [double bands] 1 bodyband-type A model & 1 B model

Rental information

You can try Ultra-Ma for 2 weeks or more ! 2 weeks: from ¥ 22,000 (tax included)~ If you purchase the same device after your rental, it can be discounted by as much as 2 months' rental fee.

For more details, see Japan Selfcare Institute®

Japan Selfcare Institute Q Search

Ultra-Ma Ultrasonic device for the head and body



This is a set of 1 headband-type ultrasonic device for the head and 1 bodyband-type ultrasonic device.

- The headband-type ultrasonic device for the head has 4 of the ultrasonic transducers used in the UM-123X model.
- The ultrasonic device for the body stimulates with ultrasonic transducers with about 10 times the intensity of the ultrasonic device for the head.
- The ultrasonic device for the body sends strong and weak waves through the body about 6,000 times/second, transmitting micro movements to your anatomy.
- Use the body band by wrapping it directly against the skin, or placing it against the area that concerns you.
- The ultrasonic device for the body comes with the A model, but if you prefer it can be changed to the B model.
- Use the ultrasonic device for the body and ultrasonic device for the head at the same time.
- You can adjust the strength of the ultrasonic waves and the volume of the music.
- It contains 3 tunes by Mozart.
- You can connect to the external port to listen to any music you like.

Set contents

Bp-21

Control unit



Headband-type

Navy blue

Ultrasonic device Ultrasonic device Music oscillator **Vibration**



〈 A model 〉

the A model or B model



Bp-21 ① ··· ¥2,200,000 (tax included) 〈 Contents 〉 1 headband-type ultrasonic device & 1 bodyband-type A model Bp-21 ② ··· ¥2,300,000 (tax included) 〈 Contents 〉 1 headband-type ultrasonic device & 1 bodyband type B model

⟨ B model ⟩

Ultra-Ma Ultrasonic device for the jaw joint



This is headphone-type ultrasonic device. You can wear it over the jaw joint and masticatory muscles.

- The headphone-type ultrasonic device stimulates with ultrasonic transducers with about 10 times the strength of the ultrasonic device for the head.
- It sends strong and weak waves through the body about 6,000 times/second.
- The length of the headphones' right and left arms is adjustable.
- You can adjust the strength of the ultrasonic waves and the volume of the music.
- It contains 3 tunes by Mozart.
- You can connect to the external port to listen to any music you like.



Price Hp-20 ··· ¥1,980,000 (tax included) ⟨ Contents ⟩ 1 headphone-type ultrasonic device



Control unit



Headphone type ultrasonic device



* The bodyband is not worn over clothing,

but rather directly against the skin.

Bodyband

* As of January 2023, prices currently reflect the 10% consumption tax, but in the event of a change in the tax, prices will also change.